IFU\_Roll\_on\_cuffs\_Ver08/2022-05-24

# Instruction

## Hpm™ roll-on cuff system

*For better care!*

SVENSKA • ENGLISH • ESPAÑOL • DEUTSCH • FRANÇAIS  
SUOMEN • ITALIANO • NORSKA • NEDERLANDS

**Symbol definition:**

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## Index

1.	Introduction	5
2.	Content	5
3.	Safety instructions	5
4.	Getting started	5
5.	Cleaning	6
6.	Sterilization	7
7.	Storage	7
8.	Liability & Warranty	7

### 1. Introduction

The Hpm™ roll-on cuff is a complete system for creating and maintaining a bloodless field in the arms and legs when performing operation in bloodless fields.

Before using the roll-on cuff system, make sure that you have received every item ordered and read through the operation instruction.

### 2. Content

The Hpm™ roll-on cuff system consists of:

- 6 roll-on cuffs in different sizes
- Hand pump with manometer
- Vacuum pump
- Hose with needle
- Wedge in rubber
- Measuring Tape

### 3. Safety instructions

The Hpm™ roll-on cuffs are delivered non sterile and should be sterilized prior to use. How to sterilize and what to consider can be read under the chapter sterilization.

**Note!** Do not sterilize hand pump and vacuum pump.

After sterilization the product must be cooled off before use. The material retains heat so even if the surface feels cold please ensure that the product is sufficiently cooled off. This is extremely important during bloodless field procedures as the skin is not cooled naturally by the blood circulation. We recommend that cuffs and rubber wedges are exposed to room temperature (20 °C) for at least 4 hours before use. Always make a safety check before use.

The roll-on cuff should not be exposed to strong solvents, sharp objects or strong UV-light as this may affect the material characteristics and product function.


**Note!** If the patient is diabetic and has been treated with cortisone be sure that the patient's skin is in good condition.

### 4. Getting started


1. Measure the circumference with the measuring tape on the patient's extremity where the roll-on cuff is planned to be placed. Choose an adequate size of roll-on cuff. We recommend using the smaller size in borderline cases.

2. Connect the hand pump and hose, please note that the hose should be sterile while the hand pump should not be sterilized.

3. Insert the needle as shown in the picture. Holding the valve according to the picture ensures that the needle does not mistakenly puncture the rubber wall once the valve has been penetrated. The valve has a concealed hole which goes straight into the cuff, perpendicular to the surface of the roll-on cuff. It is very important to insert the needle in this direction and to carefully follow the concealed hole. Following the concealed hole with the needle makes it relatively simple to insert the needle. If inserting the needle requires a lot of strength, it is a sign that you have missed the concealed hole. The needle should then be withdrawn so that the user can start finding the concealed hole again. Repeatedly inserting the needle outside the concealed hole risks destroying the valve's function and can lead to the valve leaking air.

*The needle must be inserted at right angles to the surface. Hold the valve as in the picture to avoid the needle protruding through the cuff on the back.*

*If the needle is inserted obliquely, it can lead to the cuff not remaining sealed.*

4. Pump the pressure in the roll-on cuff up to 130 mmHg.
5. Remove the needle
6. Roll the roll-on cuff up the patient's extremity. When rolling it up distally the extremity is exsanguated. The blood flow stops where the roll-on cuff is placed.
7. To ensure that the roll-on cuff stays in the desired position, a wedge may be used to secure the cuff.
8. The roll-on cuff exerts a pressure of 200 to 320 mmHg against the extremity. Due to the unique design of the roll-on cuff, the pressure will never exceed 350 mmHg. This is an integrated safety feature of the Hpm™ roll-on cuff system.
9. Carefully remove the roll-on cuff by rolling it off the extremity when the procedure and the bandaging has been completed.
10. Clean soiled cuffs and wedges by washing them with approved and tested detergent for silicon products. For satisfactory aseptic cleaning, the roll-on cuffs should be cleaned by being washed in a disinfectant for 90°C for at least 10 minutes. Do not clean the roll-on cuffs with mineral oil-based detergent. Further information about this is available in chapter 5.
11. Sterilise the roll-on cuffs in an autoclave at 134°C for 3 minutes. Further information about this is available in chapter 6.

## 5. Cleaning

To clean the roll-on cuff use solvents approved for silicone materials. Hard brushes and other cleaning methods which can cause damage to the surface should never be used. Any cleaning agent residue should be flushed off carefully with running water.

The roll-on cuff can be washed in a dishwasher at a maximum highest temperature of 90 °C.

### **Do not clean the roll-on cuffs with mineral oil based cleaning agents.**

Avoid using cleansing agents that contain ozone, because the ozone breaks down the material much faster than the products normal lifetime.

**We recommend cleaning in disinfectant rather than hand washing.**

## 6. Sterilization

1. Empty the roll-on cuff by deflating it with the vacuum pump.  
It is essential to empty all the air from the roll-on cuff.

**The roll-on cuff might burst in the autoclave during the vacuum phase if there is air left in it. If the roll-on cuff should burst it contains 1ml of silicon oil.**

2. Autoclave at 134 °C for 3 minutes.  
**We recommend sterilizing it with steam.**  
**Warning!** Do not sterilize the roll-on cuff with Ethylene Oxide. Do not sterilize the roll-on cuff in dry heat sterilizers.

3. The hose and needle are sterilized in the same way as the roll-on cuffs.  
The measuring tapes must be sterilized with the surfaces kept apart (i.e. with gauze bandage) to avoid loss of the color.

**Very Important:** The roll-on cuffs and wedges must be properly cooled off before they are used on patient.

The material used in the roll-on cuff system retains heat so even if the surface feels safe to touch it is necessary to ensure that all heat has been transmitted away before use. This is especially important in the bloodless field because in this situation the skin is not cooled by blood circulation. **Therefore, it is recommended that sterilization is carried out at least 4 hours before use, allowing the components to cool at room temperature (20 °C).** Do not attempt to reduce the cooling off time by e.g. dipping the wedge in sterile water. A final check must always be made immediately prior to use to be absolutely certain that the wedges and roll-on cuffs have been completely cooled.

The product is designed to be used at least 25 times over a maximum period of 2 years

## 7. Storage

We recommend that the product is packed in a light-proof plastic bag. The products must not be exposed to UV-light and temperatures outside the range 10-25 °C.

Hpm™ roll-on cuffs are delivered in a black plastic bag in order to prevent exposure to UV-light. The roll-on cuff will have a longer lifetime if stored in a sealed black plastic bag. Sterilization in autoclave with steam does not deteriorate the silicone as much as UV-light and ozone. We recommend that the product is stored in a dark environment.

## 8. Liability & Warranty

In order to prevent product malfunction, we recommend replacing the individual roll-on cuffs according to the warranty information.

The information in this instruction is not compulsive and does not deprive the user from the responsibility of checking the components of the system before use. The manufacturing warranty does not cover damage caused through misuse or normal wear and tear.

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Hammarplast Medical AB