

Hpm™

Instruction Tablet crusher Manual Tablettkross

Artno. 300280, 300320

EN Cleaning

Dishwasher
max 70° C or
wipe with a
cloth and
disinfectant
solution.

NOTE! The
crusher can
not be auto-
claved!

Warranty

1 year from
date of
purchase.

SE Rengöring

I diskmaskin
max 70° C alt.
torka av med
trasa och
desinfektions-
sprit.


OBS! Krossen tål
ej att autokla-
veras!

Garanti


1 år från
inköpsdatum.




1




2



3



4



Don't forget to
clean the crusher.



Hpm[™]

Hammarplast Medical AB
Kartåsgatan 8, SE-531 40 Lidköping, Sweden
Tel. +46(0)510 61880, info@hpm.se
www.hpm.se