

ECO+ Disposable Fossil Free Ear Speculum



Product information

Ref. No.	Name	Colour	Diameter	Height
31210	ECO+ ear speculum ø 2,5 mm	Black	ø 2,5 mm	36,4 mm
31220	ECO+ ear speculum ø 4 mm	Black	ø 4 mm	34 mm

Packaging

Ref. No.	Inner packaging	Outer packa- ging	Car- tons/ pallet	Carton size mm	Pallet height	Carton weight	Pallet weight	GTIN carton	GTIN sleeve
31210	250 pcs/ box	7500	32	585x360x200	180 cm	4 kg	152 kg	7340077904096	7340077905024
31220	250 pcs/ box	7500	32	585x360x200	180 cm	4 kg	158 kg	7340077904102	7340077905031

Information

Hpm™ disposable ear speculums are accessories for otoscopes used for direct examination of the external auditory canal and the tympanic membrane.

- Available in two different sizes: ø 2,5 mm och ø 4 mm.
- Made from fossil-free material (PP).
- Disposable ear speculums are hygienic, safe and prevent cross-contamination.
- Secure fit.
- No sharp edges.
- Available in black.
- Dp not contain BPA, phthalates, latex or PVC.
- Hpm™ ear speculums are CE-marked and classified as a Class 1 medical device.
- For single-use.

Labels

ECO+

EN Disposable Ear Speculum
DE Einmal Ohrtrichter
FR Spéculums auriculaires jetables
SE Örontrattar för engångsbruk
ES Espéculo auriculares desechable
DK Øretragt engangs
NO Øretube engangs
FI Kertakäyttöiset korvasuppilot

ø 2,5 mm

Svart/Schwarz/Noir/Svart/Negro/Sort/Svart/Musta

Pcs/Stck/Pcs/St/Stk/Kpl: 30x250

REF 31210


Hammarplast Medical AB
Kartåsgatan 8
531 40 Lidköping
Sweden
info@hpm.se

   2028-09

  

 7 340077 904096

Ver nr 9/2023-04

    2028-09

   230905
15:48

 **Hammarplast Medical AB**
Kartåsgatan 8, 531 40 Lidköping, Sweden
info@hpm.se

  **ISCC**
International Sustainability
& Carbon Certification

REF 31210

 7 340077 904096

Ver nr 10/2023-04

Symbols



Manufacturer



Do not re-use



Product conforms with applicable EC directives. Notified Body-number. Inspected yearly by Intertek.



Keep dry



Reference number



Made from renewable raw material. Inspected yearly by the independent organisation ISCC.



Keep away from sunlight



Batch/ LOT-number



Medical Device



Use-by date



Medical Device

Hammarplast Medical AB
Kartåsgatan 8
SE-531 40 Lidköping, Sweden

Tel +46 (0)510-618 80
www.hpm.se

HpmTM

For better care!