

# Gallipot 250 ml



## Product information

Ref. No.	Product name	Color	Diameter	Height	Volume
200400	Gallipot	Transparent	94,40 mm	48,55 mm	250 ml
200410	ECO+ Gallipot fossil free	Transparent	94,40 mm	48,55 mm	250 ml

## Packaging

Ref. No.	Inner packaging (sleeve)	Outer packaging (carton)	Cartons/pallet	Carton size	Pallet height	Carton weight	Pallet weight	GTIN carton
200400	16 pcs	240 pcs	60 pcs	398x198x368 mm	200 cm	2,7 kg	182 kg	7340077901408
200410	16 pcs	240 pcs	60 pcs	398x198x368 mm	200 cm	2,7 kg	182 kg	7340077901422

## Information


**Gallipot is a disposable product with a capacity of 250 ml. It is useful as a container for sampling, short-term storage, microanalysis, and many other applications.**


- Temperature resistant up to 120 °C
- Transparent with a clear scale
- Volume 250 ml
- Can be sterilized
- Available in PP and fossil-free material

# Labels

**Gallipot 250 ml**  
Transparent  
**Gallipot 250 ml**  
Transparent  
Pcs/Antal: 240 (15 pcs x16 )  
**ART. NR: 200400**


**Hpm**<sup>TM</sup>  
*For better care!*

  
**Hammarplast  
Medical AB**  
Kartåsgatan 8, SE-531 40  
Lidköping, Sweden  
www.hpm.se, info@hpm.se

  
7 340077 901408  
Vers. 01/2023-03



**LOT** 250225 14:48:52


EN0002


  
**ECO+ Gallipot 250 ml**  
Transparent  
**ECO+ Gallipot 250 ml**  
Transparent  
Pcs/Antal: 240 (15 pcs x16 )  
**REF 200410**


**Hpm**<sup>TM</sup>  
*For better care!*

**LOT** 250225

  
International Sustainability  
& Carbon Certification

  
7 340077 901422  
Vers. 01/2023-03

  
**Hammarplast  
Medical AB**  
Kartåsgatan 8, SE-531 02  
Lidköping, Sweden  
www.hpm.se, info@hpm.se

EN0002

# Symbols



Manufacturer



Reference  
number



Keep dry



Batch/LOT-number



Keep away from  
sunlight



Made of renewable raw  
material, inspected yearly  
by the independant organi-  
sation ISCC.

**Hammarplast Medical AB**  
Kartåsgatan 8  
SE-531 40 Lidköping, Sweden

Tel +46 (0)510-618 80  
[www.hpm.se](http://www.hpm.se)

**Hpm**<sup>TM</sup>

*For better care!*