

- ISO 13485
- ISO 9001
- ISO 14001

Pill Straw





Product information

Ref. No.	Description	Colour	Top diameter	Bottom diameter	Height
300700	ECO+ Pill straw in carton 2-p	Blue		Outer: 9 mm Inner: 7 mm	110 mm
300701	ECO+ Pill straw in bag 10-p	Blue	Outer: 14 mm Inner: 11 mm	Outer: 9 mm Inner: mm	110 mm

Packaging

Ref. No.	Inner package	Outer package	Outer cases per pallet	Inner case size mm	Outer case size mm	Package height	Pallet height	Weight package	Weight pallet	GTIN
300700	2 pcs	20 pcs	240	160x50x20	180x130x120	120 cm	1 040 cm	0,2 kg	63 kg	7340077902115
300701	10 pcs	80 pcs	240	200x120x30	180x130x120	120 cm	1040 cm	0,3 kg	87 kg	7340077902122


Information

Pill Straw makes intake of tablets easier.

- Place the tablet inside the straw and position the straw on the glass.
- A clip keeps the straw securely in place.
- Can be reused multiple times.
- Made from fossil-free materials.

To ensure a high standard of hygiene and ease of use, we recommend the following cleaning instructions:

- After use, make sure to clean the Pill Straw to remove any residue.
- The Pill Straw is dishwasher safe at 60°C.
- For best results, use standard dishwashing detergent and disinfectant.
- After machine washing, rinse the straw thoroughly with clean water to ensue no detergent residues remain inside.



Labels





Pill Straw 2-p

Pill Straw 2-p

Blue

Pcs/Antal: 10x2 pcs

REF 300700

Hammarplast Medical AB

Kartåsgatan 8, SE-531 40 Lidköping, Sweden www.hpm.se, info@hpm.se



Vers. 01/2023-09

LOT 240821







Pill Straw 10-p

Pill Straw 10-p

Blue

Pcs/Antal: 8x10 pcs

REF 300701

Hammarplast

Kartåsgatan 8, SE-531 40 Lidköping, Sweden www.hpm.se, info@hpm.se



LOT 240821

Pill Straw 10-p

LOT

240821 13:59

Hammarplast Medical AB

Kartåsgatan 8. 531 40 Lidköping, Sweden info@hpm.se



REF

300701



Vers:1/2023-12

Symbols



Manufacturer



Made of renewable raw material, inspected yearly by the independant organisation ISCC.



Reference number



Batch/LOT-number

Hammarplast Medical AB Kartåsgatan 8 SE-531 40 Lidköping, Sweden

Tel +46 (0)510-618 80 www.hpm.se

